# BC 58

(D) **Blutdruckmessgerät**
Gebrauchsanweisung ..... 2–14
(GB) **Blood pressure monitor**
Instructions for use ..... 15–26
(F) **Tensiomètre**
Mode d'emploi ..... 27–39
(E) **Tensiómetro**
Manual de instrucciones ..... 40–52
(I) **Misuratore di pressione**
Istruzioni per l'uso ..... 53–65
(TR) **Bilgisayarli tansiyon ölçer**
Kullanım kılavuzu ..... 66–77
(RUS) **Прибор для измерения давления на запястье**
Инструкция по применению ..... 78–92
(PL) **Ciśnieniomierz obsługi**
Instrukcja obsługi ..... 93–105
**Electromagnetic Compatibility Information** ..... 106–108

BEURER GmbH • Söflinger Str. 218 • 89077 Ulm (Germany)
Tel.: +49 (0) 731 / 39 89-144 • Fax: +49 (0) 731 / 39 89-255
www.beurer.com • Mail: kd@beurer.de

CE 0483

# DEUTSCH

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# ENGLISH

1. Getting to know your instrument ..... 15
2. Important information ..... 16
3. Unit description ..... 19
4. Prepare measurement ..... 20
5. Measuring blood pressure ..... 21
6. Evaluating results ..... 22
7. Saving, retrieving and deleting results ..... 23
8. Transferring measurements ..... 24
9. Error messages/trouble-shooting ..... 24
10. Cleaning and storing the instrument ..... 25
11. Specifications ..... 25

**Dear Customer,**

thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team



## 1. Getting to know your instrument

The wrist blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure.
You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings together with the average.
A warning is issued for anyone suffering from cardiac arrhythmia.
The values determined are classified and graphically evaluated according to WHO guidelines.
Keep these instructions for use carefully for further use and also let other users have access to them.

## 2. Important information

### Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

| | |
|---|---|
|  | Caution |
|  | Note<br>Note on important information |
|  | Follow instructions for use |
|  | Type BF applied part |
|  | Direct current |
|  | Please dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). |
|  | Manufacturer |
|  | Permissible transport and storage temperature. Permissible transport and storage humidity. |
|  | Permissible operating temperature and humidity. |
|  | Keep dry |
| SN | Serial number |
| CE 0483 | The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products. |

### Advice on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes after eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.

- The measurements taken by you are for your information only – they are not a substitute for a medical examination! Discuss the measurements with your doctor, and never base any medical decisions on them (e.g. medicines and their administration)!
- Do not use the blood pressure monitor on newborns, pregnant women or patients with preeclampsia.
- In the case of restricted circulation on the arm as a result of chronic or acute vascular diseases (including vascular constriction), the accuracy of the wrist measurement is limited. In this case you should avoid using an upper arm blood pressure monitor.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the unit on people who have the specified wrist measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- The blood pressure monitor can only be operated with batteries.
- To conserve the batteries, the monitor switches off automatically if no buttons are pressed for one minute.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

### ⚠ Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
  – You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.

– Never drop the device.
– Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
– Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.

- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

### ⚠ Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, flush out the affected areas with water and seek medical assistance.
- ⚠ **Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.
- ⚠ **Risk of explosion!** Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only .
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.

### ⓘ Repair and disposal

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Beurer Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

## 3. Unit description


1. Display
2. Memory button **MEM**
3. START/STOP button ⏻
4. USB interface
5. Touch screen activation switch
6. Battery cover

**Touch screen activation switch**

The device features a touch screen display. To prevent the screen from being activated accidentally, keep the touch screen activation switch set to **OFF** when the device is not in use. To operate the device, move the touch screen activation switch to **ON**. When the touch screen display is touched (START/STOP button ⏻ or **MEM** button), a beep is output.

**Tip:** You can switch off the device at any time by setting the touch screen activation switch to the **OFF** position.

**Icons in the display:**

1. Time and date
2. WHO classification
3. Systolic pressure
4. Diastolic pressure
5. Measured pulse
6. „Change battery“ icon
7. Pump up, release air (arrow)
8. Memory display (M), average value (A), morning (AM), evening (PM)
9. Cardiac arrhythmia icon ((♥))

**PC interface**

The Beurer blood pressure monitor also allows you to transfer your measured values to the PC.
To do this you require a USB cable (included in delivery) and the Beurer “Health Manager” PC software.
The software can be downloaded free of charge from www.beurer.com/service/download

**System requirements for the Beurer “Health Manager” computer software**

**1. Supported operating systems:**


- Windows XP SP3
- Windows Vista SP1 or later
- Windows 7
- Windows 7 SP1

**2. Supported architectures:**
- x86 (32 bit)
- x64 (64 bit)

**3. Hardware requirements:**
- Recommended: at least Pentium 1 GHz or faster with at least 1 GB RAM
- Free memory on the primary partition of at least:
  – x86 – 600 MB
  – x64 – 1.5 GB
- Graphic resolution from: 1024 x 768 pixels
- USB port 1.0 or later

## 4. Prepare measurement

**Inserting battery**

- Open the battery compartment lid.
- Use exclusively brand-name batteries of the type: 2 x 1.5V Micro (Alkaline Type LR 03). Making absolutely sure that you insert them with the correct polarity as marked. Do not use rechargeable batteries.
- Replace the battery cover carefully.

If the battery change is continuously illuminated, measurement is no longer possible and you must replace all the batteries. Once batteries have been removed from the device, the time must be reset.

The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.
The codes below are printed on batteries containing harmful substances:
Pb = Battery contains lead,
Cd = Battery contains cadmium,
Hg = Battery contains mercury.

**Setting date and time**

You should set the date and time without fail. Otherwise, you will not be able to save your measured values correctly with a date and time and to access them again later.
**Tip:** If you press and hed the **MEM** button, you can set the values more quickly.


To set the date and time, proceed as follows:
- Move the touch screen activation switch to the **ON** position.
- Press the START/STOP and **MEM** buttons simultaneously, 24h begins to flash. Select 12h or 24h mode using the **MEM** button. Press START/STOP to confirm your selection. The year display will start to flash. Set the year with the **MEM** button and confirm with START/STOP ⏻.
- Then set the month, day, hour and minutes and confirm each setting with START/STOP ⏻.

- Pressing the START/STOP button ⏻ again will switch the display off.

## 5. Measuring blood pressure

Please ensure the device is at room temperature before measuring.

**Positioning cuff**

- Bare your left wrist, making sure that the circulation in the arm is not restricted by any clothes etc. that are too tight. Position the cuff on the inside of your wrist.
- Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approx. 1 cm below the ball of your thumb.
- The cuff has to be fitted tightly around the wristbut should not constrict it.

**Important:** The instrument should only be operated with the original cuff.

**Correct posture**

- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement while sitting or lying. To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground. Make sure to rest your arm and move it. Always make sure that the cuff is at heart level. Otherwise significant deviations can occur. Relax your arm and the palm of your hand.

- In order not to distort the result, it is important to keep still during the measurement and not talk.

**Measuring blood pressure**

- Move the touch screen activation switch to the **ON** position.
- Switch on the blood pressure monitor with the START/STOP button ⏻. After the full-screen display, the most recently used user memory appears (U1 or U2). To change the user memory, press the **MEM** button and confirm your selection with the START/STOP button ⏻. If no button is pressed, the most recently used user memory is automatically used after 5 seconds.

- Before the measurement, the last saved test result is briefly displayed. If there is no measurement in the memory, the instrument always displays the value 0.
- The cuff automatically inflates. Cuff air pressure is released slowly. If a tendency towards high blood pressure is already detectable, the cuff is pumped up again and cuff pressure increased further. As soon as a heart rate is detected, the heart rate symbol ♥ is displayed.
- After the pressure has been completely reduced, the systolic pressure, diastolic pressure and pulse readings are displayed.
- You can cancel the measurement at any time by pressing the START/STOP button ⏻ or by pushing the touch screen activation button to **OFF**.
- Er appears if it has not been possible to perform the measurement properly. Observe the section in these instructions for use on error messages/troubleshooting and repeat the measurement.
- The test result is saved automatically.
- To switch off the device, press the START/STOP button ⏻ or push the touch screen activation switch to **OFF**. If you forget to switch off the device, it switches off automatically after approx. 1 minute.

Wait at least 5 minutes before taking another measurement!

## 6. Evaluating results

**Cardiac arrhythmia:**

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon ♥.
This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor.
Repeat the measurement if the flashing icon ♥ is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon ♥ appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

**WHO classification:**

In accordance with the guidelines/definitions of the World Health Organization and the latest findings, the measurements can be classified and assessed according to the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

| Range of blood pressure values | Systolic (in mmHg) | Diastolic (in mmHg) | Measure |
|---|---|---|---|
| Grade 3: Severe hypertension | ≥180 | ≥110 | Seek medical advice |
| Grade 2: Moderate hypertension | 160–179 | 100–109 | Seek medical advice |
| Grade 1: Mild hypertension | 140–159 | 90–99 | Have it checked regularly by doctor |
| High-normal | 130–139 | 85–89 | Have it checked regularly by doctor |
| Normal | 120–129 | 80–84 | Check it yourself |
| Optimal | <120 | <80 | Check it yourself |

Source: WHO, 1999

The WHO classification in the display shows the range of the measured blood pressure.
If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the grade hypertension range and diastolic pressure in the normal range), the WHO classification on the unit always indicates the higher range (grade hypertension in the example described).

## 7. Saving, retrieving and deleting results

- The results of each successful measurement are stored together with date and time. With more than 60 items of measured data, the earliest items of data measured are lost.
- Move the touch screen activation switch to the **ON** position.
- Make a selection with the **MEM** button, and then confirm the desired user memory with the START/STOP button ⏻. If you press the **MEM** button again, the average value Ø of all the stored measured values in the user memory will be displayed. If you press the **MEM** button again, the average value of the last 7 days for the morning measurement will be displayed (Morning: 5 a.m. – 9 a.m., display Ø☀). If you press the **MEM** button again, the average value of the last 7 days for the evening measurement will be displayed (Evening: 6 p.m. – 8 p.m., display Ø☾). If you continue to press the **MEM** button, the most recent

individual measured values with date and time are displayed in turn.
- You can delete the memory by pressing and holding the **MEM** button for 3 seconds. All the values in the current user memory are deleted after three beeps are output.
- To switch off the device, press the **MEM** button again or the START/STOP button ⏻ or push the touch screen activation switch to **OFF**.
- If you forget to switch off the device, it will switch off automatically after 2 minutes.

## 8. Transferring measurements

Connect the blood pressure monitor to your PC using the USB cable.

ⓘ No data transfer may be launched whilst performing a measurement.

PC is shown on the display. Begin the data transfer in the "HealthManager" PC software. During the data transfer, an animation is shown on the display. A successful data transfer is displayed as in figure 1.

figure 1

If the data transfer is unsuccessful, an error message appears as in figure 2. In this case, interrupt the PC connection and start the data transfer again.

figure 2

After 30 seconds of not being in use or if communication with the PC is interrupted, the blood pressure monitor switches itself off automatically.

## 9. Error messages/trouble-shooting

In case of faults, the Er_ message appears in the display. Error messages can occur when
- systolic or diastolic pressure could not be measured (Er 1 or Er2 appears on the display),
- systolic or diastolic pressure was outside the measurement range (Hi or Lo appears on the display),
- the cuff is fastened too tightly or loosely (Er3 or Er4 appears on the display),
- the pump pressure is higher than 300 mmHg (Er5 appears on the display),
- pumping up takes longer than 160 seconds (Er6 appears on the display),
- there is a system or device error (ErA, Er0, Er7 or Er8 appears on the display),
- the batteries are almost empty ,
- the data could not be sent to the PC (PC Er appears in the display).

In such cases, repeat the measurement. Ensure you do not move or speak.
If necessary, re-insert or replace the batteries.

**Technical alarm – description**

Should the recorded blood pressure (systolic or diastolic) lie outside the limits specified in the section “Technical specifications”, the technical alarm will appear on the display indicating either “Hi” or “Lo”. In such cases, you should seek medical assistance and check the accuracy of your procedure.

The limit values for the technical alarm are factory set and cannot be adjusted or deactivated. These alarm limit values are accorded second priority under the standard IEC 60601-1-8.

The technical alarm is a non-locking alarm and must not be reset. The signal shown on the display will disappear automatically after about 8 seconds.

## 10. Cleaning and storing the instrument

- Clean your blood pressure monitor carefully using a slightly damp cloth only.
- Do not use detergents or solvents.
- On no account must you immerse the computer in water, otherwise liquid can enter it and cause demage.
- When storing the device, make sure that no heavy objects are placed on top of it. Remove the batteries. The cuff tube should not have any sharp kinks.

## 11. Specifications

| | |
|---|---|
| Model no. | BC 58 |
| Measurement method | Oscillometric, non-invasive blood pressure measurement on the wrist |
| Measurement range | Cuff pressure 0–300 mmHg,<br>systolic 60–260 mmHg,<br>diastolic 40–199 mmHg,<br>Pulse 40–180 beats/minute |
| Display accuracy | Systolic ±3 mmHg, diastolic ±3 mmHg, pulse ±5 % of the value shown |
| Measurement inaccuracy | Max. permissible standard deviation according to clinical testing:<br>systolic 8 mmHg/diastolic 8 mmHg |
| Memory | 2 x 60 memory spaces |
| Dimensions | L 90 mm x W 68 mm x H 30 mm |
| Weight | Approx. 158 g (without batteries) |
| Cuff size | 140 to 195 mm |
| Permissible operating conditions | +5 °C to +40 °C, ≤90 % relative air humidity (non-condensing) |
| Permissible storage conditions | -20 °C to +55 °C, ≤95 % relative air humidity, 800–1050 hPa ambient pressure |
| Power supply | 2 x 1,5V ⎓ AAA batteries |

| | |
|---|---|
| Battery life | For approx. 300 measurements, depending on the blood pressure level and/or pump pressure |
| Accessories | Instruction for use, 2 x 1.5 V AAA batteries, storage box |
| Classification | Internal supply, IPX0, no AP or APG, continuous operation, type BF applied part |

Technical information is subject to change without notification to allow for updates.

- This unit is in line with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Service address or found at the end of the instructions for use.
- This device is in line with the EU Medical Devices Directive 93/42/EC, the “Medizinproduktegesetz” (German Medical Devices Act) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2 – 30: Particular requirements for the safety and essential performance of automated non-invasive blood pressure monitors).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

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# Electromagnetic Compatibility Information

**Table 1**
**For all ME EQUIPMENT and ME SYSTEMS**

| *Guidance and manufacture's declaration – electromagnetic emissions* | | |
|---|---|---|
| The BC 58 is intended for use in the electromagnetic environment specified below.<br>The customer or the user of the BC 58 should assure that it is used in such an environment. | | |
| **Emissions test** | **Compliance** | **Electromagnetic environment – guidance** |
| RF emissions CISPR 11 | Group 1 | The BC 58 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The BC 58 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Not applicable | |

**Table 2**
**For all ME EQUIPMENT and ME SYSTEMS**

| *Guidance and manufacturer's declaration – electromagnetic immunity* | | | |
|---|---|---|---|
| The BC 58 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 58 should assure that it is used in such an environment. | | | |
| **IMMUNITY test** | **IEC 60601 test level** | **Compliance level** | **Electromagnetic environment – guidance** |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact<br>±8 kV air | ±6 kV contact<br>±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| **NOTE:** $U_T$ is the a.c. mains voltage prior to application of the test level. | | | |

**Table 3**
**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

*Guidance and manufacturer's declaration – electromagnetic immunity*

The BC 58 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 58 should assure that it is used in such an environment.

| IMMUNITY test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
|---|---|---|---|
| | | | Portable and mobile RF communications equipment should be used no closer to any part of the BC 58, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. |
| | | | **Recommended separation distance:** |
| Conducted RF IEC 61000-4-6 | 3Vrms 150 kHz to 80 MHz | 3V | $d = 1.2 \sqrt{P}$ |
| | | | $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz |
| Radiated RF IEC 61000-4-3 | 3V/m 80 MHz to 2.5 GHz | 3V/m | $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz |
| | | | Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m). |
| | | | Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,[a] should be less than the compliance level in each frequency range.[b] Interference may occur in the vicinity of equipment marked with the following symbol:  |

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.
**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

[a] Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BC 58 is used exceeds the applicable RF compliance level above, the BC 58 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BC 58.

[b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

**Table 4**

**For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

*Recommended separation distances between portable and mobile RF communications equipment and the BC 58*

The BC 58 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BC 58 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BC 58 as recommended below, according to the maximum output power of the communications equipment.

| **Rated maximum output power of transmitter** (W) | **Separation distance according to frequency of transmitter** (m) | | |
|---|---|---|---|
| | **150 kHz to 80 MHz** $d = 1.2\sqrt{P}$ | **80 MHz to 800 MHz** $d = 1.2\sqrt{P}$ | **800 MHz to 2.5 GHz** $d = 2.3\sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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